Super Rapid Readout Fluorescence System

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Parameters determined at time of manufacture according to ISO 11138-1: 2006 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

ISO and USP Compliant
 ATCC is a registered trademark of American Type Culture Collection

Lic. Adrián J. Rovetto
 Director Técnico
 Technical Director

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Biological Indicators

For plasma or vaporized Hydrogen Peroxide sterilization

Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC 7953 spores inoculated on a strip (spores carrier). It also has a purple colored culture media contained within the glass ampoule, next to the carrier, on the base of the tube.

Product description

Bionova® BT96 Super Rapid Readout Biological Indicators have been designed for quick and easy monitoring of Plasma or Vaporized Hydrogen Peroxide sterilization cycles. They consist of a plastic tube, a spore carrier and a glass ampoule with a culture medium. The tube has a plastic top with holes and a barrier permeable to Hydrogen Peroxide.

Super Rapid Readout: 30 minutes

The Super Rapid Readout must be carried out in the Bionova® IC10/20FR, IC10/20FRCLD or MiniBio Reader Incubator or in a similar Reader Incubator that meets the features that are described below. The Reader Incubator must be capable of reading the fluorescence emission of the product, resulting from breaking a specific substrate on the spores carrier. Fluorescence is produced when the reader excites the carrier with UV light at 360 nm. Final readout of negative results is readily available at 30 minutes of incubation. Fluorescence is a direct process of the germination and growth of *Geobacillus stearothermophilus* spores which have survived the sterilization process (positive result). Furthermore, a failure in the sterilization process can also become evident by the color change of the culture media. Due to the high sensitivity of the fluorescence results after 30 minutes, conventional incubation for color change of Bionova® BT96 Biological Indicator is not an advantage.

ATTENTION: Do not use the Bionova® BT96 Biological Indicators to control sterilization process for Ethylene Oxide, Dry Heat, Formaldehyde, Steam or other sterilization processes different from the specified in this IFU. Do not reuse the biological indicators.

Storage

Best stored in the original box under the following conditions: temperature between 10-30 °C, relative humidity between 30-80%.

Do not freeze.

Do not store biological indicators near sterilizing agents or other chemical products.

Instructions for use

1. Identify the Bionova® BT96 Biological Indicator by writing the sterilizer number (in case of having more than one), load number and processing date on the label.
2. Pack the biological indicator along with materials to be sterilized in an appropriate package according to recommended sterilization practices. Place the package in those areas that you consider a priori most inaccessible for the sterilizing agent (hydrogen Peroxide). Generally, the center of the load and areas near the doors are more problematic.
3. Sterilize as usual.
4. After the sterilization process has finished, remove the biological indicator from the package.
5. Check that the chemical indicator printed on the label of the biological indicator changes to green. A color change indicates that the biological indicator has been exposed to Hydrogen Peroxide. **IMPORTANT:** This color change does not evidence the process effectiveness to achieve sterility. If the chemical indicator color has not changed, check the sterilization process.
6. Start the Bionova® IC10/20FR, IC10/20FRCLD or MiniBio Reader Incubator, select the proper temperature (60°C) and incubation program (30 minutes). **IMPORTANT:** Do not reveal the Bionova® BT96 Biological Indicator until the incubator reaches the required temperature.
7. Press the lid to seal the tube. Then, break the ampoule contained in the biological indicator. This may be done in 3 different ways:
 A- With an ampoule crusher.
 B - Manually. To do this, take the tube by placing your index finger and thumb on the space between the cap and the line of chemical indicator and then compress the tube.
 C- With the ampoule crusher placed within the top of the incubator's incubation area. Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the media reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator.

IMPORTANT: Use a non-sterilized biological indicator as a positive control every time a processed indicator is incubated. The positive control helps ensure that correct incubation conditions were met; capability of the media to promote rapid growth; viability of spores has not been altered due to improper storage temperature, humidity or proximity to chemicals and proper functioning of Bionova® IC10/20FR, IC10/20FRCLD o MiniBio Reader Incubator. Both the positive control and the processed indicator should belong to the same batch.

8. Incubate the processed biological indicator and the indicator used as positive control for a maximum of 30 minutes at 60±2 °C for Super Rapid Readout. Fluorescence detection by the reader (excitation 340-380 nm / emission 455-465 nm) means a failure in the sterilization process. If no fluorescence is detected at the end of the 30-minute incubation, then the result is negative (the sterilization process has been effective). The indicator used as positive control must be detected as such by the reader. It is recommendable to incubate a positive control for a visual color change. Record the positive results and discard the biological indicators immediately, as shown below.

WARNING: Do not reuse the sterilizer until the biological indicator test result is negative. Once the fluorescence readout of the biological indicators has finished, do not restart it, since results might be erroneous.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the sterilization process has not been successful, culture media will change to yellow during incubation at 60 °C, indicating the presence of living spores. If sterilization was successful, culture media will remain purple after the incubation process.

A definitive negative result is obtained after the 48-hour incubation. The positive control must show a color change from purple to yellow for results to be valid.

A definitive negative result is obtained after the 48-hour incubation. The positive control must show a color change from purple to yellow for results to be valid.

Visual readout: 7 days

A 7-day readout is optional and not intended to be routinely performed. This is an initial validation of the 30-minutes readout. The fluorescence results are compared to the 7-day visual readouts. The sensitivity of the system is determined as the difference between the 7-day positive indicators and false negative indicators (negative fluorescence reading and visually positive) in relation to the positives at 7 days (*). Based on the claimed ≥ 97% readout reliability of the 30-minutes biological indicator, there is no advantage to include the indicator beyond 30 minutes.

(* NOTE: If 7-day readout is performed, a humidified environment will be required to prevent media from drying out.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved at 121 °C for at least 20 minutes or at 132 °C for 15 minutes in a gravity displacement Steam sterilizer, or at 134 °C for 10 minutes in a vacuum assisted Steam sterilizer.