
VIEW - TEST PACK®

PAQUETE BOWIE & DICK DESECHABLE BOWIE-DICK TEST PACK

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BOWIE & DICK TEST PACK INSTRUCTIONS SHEET

0. *Introduction*

Working test for autoclaves provided with air removal system (Bowie & Dick test) has been performed to demonstrate the efficacy of non-condensable gases removal from the chamber, its leakproof capacity and the sterilizing agent capacity to penetrate the loads.

1. *Materials*

- 1.1. A sterilizer with water steam system working with air removal system on the previous phase.
- 1.2. Test Pack.

2. *Test achievement*

- 2.1. Test shall be carried out daily before beginning the sterilization routine work.
- 2.2. First of all, a cycle is run with the empty chamber to prepare it and to pick up any possible condensable gases present in the entry tubing.
- 2.3. Then, the test pack is placed inside the sterilizer as an unique load (other material has not to be placed inside the sterilizer), making sure that the upper side of the pack is upwards.
- 2.4. Place the pack on a tray at the coolest zone of the chamber, this is, down inside the autoclave, between the drain and the door.
- 2.5. Close the door and run a sterilization cycle, avoiding drying time making longer. In equipments provided with a specific cycle for Bowie & Dick, this one will be used.

3. *Results interpretation*

- 3.1. Once the sterilization process is over, remove the pack from the sterilizer and examine the colour change on the test sheet, visually.
- 3.2. Indicator colour turns to DARK BROWN-BLACK.
 - 3.2.1. If colour change is uniform, the autoclave working conditions are satisfactory. CORRECT.
 - 3.2.2. If colour change is not uniform or, by visual examination, different colour shades are observed among different zones of the indicator, a failure is being detected on the equipment working. INCORRECT.
 - 3.2.3. If colour changes to gray shades of metallic aspect, sterilization conditions exceed either time or temperature parameters. CORRECT.

4. *Storage*

- Store under room temperature (between 10 and 30°C)
- Protect from excessive humidity
- Do not use packs after its expiration date.
- Sheets neither change nor turn back its colour after sterilization process and do not need special storage conditions.
- Record relevant test cycle information, in space provided on the test sheet which should be filed with the cycle graphic record.